

DOLOMITE | JAZZ



An innovation in rollator design

DOLOMITE JAZZ combines modern design features with exceptional functionality. The clever cross folding design mechanism means it is quick and simple to fold away for storage or transportation. It offers various height and handle adjustments combined with a built-in brake system to provide stability, comfort and safety for a wide range of individuals.





Easy to operate, simple to store

Made of aluminium, the lightweight **DOLOMITE JAZZ** folds simply into a compact, upright position allowing for easy storage. The spacious basket is made of a durable and flexible material that folds with the rollator. It offers an easy to use brake system that is integrated into the frame ensuring a sleek and stylish design.

Adaptable to individual needs

The **DOLOMITE JAZZ** is available in two heights with easily adjustable handles that can be set in precisely the right position to suit individual needs. The seat has been designed for comfort and a back support comes as standard.

NEW! Slow Down Brake - A safety feature that adds an adjustable resistance to one or both wheels as required.

Features and Options



Comfortable backrest

A comfortable and safe support to lean against. Mounted in the handle bars.



Cane holder

Easy to fit, easy to use. No tools.
Part no 1512912



Basket Cover

Stylish and simple to fit, gives the basket weather protection and privacy. Comes in two colours - red or silver.
Part No 1524818 (Red)
Part No 1524819 (Silver).



Slow-down brake

A safety feature that prevents the walker rolling away. Perfect for gait training and one handed solutions.
Part no 13041-21

Technical data

JAZZ 600	620 mm	670 mm	440 mm	750 - 970 mm	600 mm
JAZZ 520	620 mm	670 mm	440 mm	670 - 810 mm	520 mm
JAZZ 600	390 mm	200 mm	7,8 kg	130 kg	260 mm
JAZZ 520	390 mm	200 mm	7,2 kg	130 kg	260 mm

Colour



Charcoal Grey

Product Code

JAZZ 600 1507823
JAZZ 520 1507775

Including tray, basket and back-strap.

INVACARE® Limited
Pencoed Technology Park - Pencoed -
Bridgend CF35 5HZ
Tel: +44 1 656 776 222
Fax: +44 1656 776 220
E-mail: uk@invacare.com
www.invacare.co.uk
Sales Order E-mail:
ordersuk@invacare.com

INVACARE® Ireland Ltd
Unit 5 - Seatown Business Campus
Seatown Road - SWORDS - County
Dublin
Tel. +353 1 8107084
Fax +353 1 8107085
Email: Ireland@invacare.com
www.invacare.ie

© 2010 Invacare® International Sarl
All rights reserved.

All information quoted is believed to be correct at time of print.
Invacare® reserves the right to alter product specifications without prior consultation.

Jazz - UK - 11/2010 - 1504398



Yes, you can.®

